



Phone: (215) 357-3911

FAX: (215) 357-5232

WEB: Abraxiskits.com

## Microcystins-LW Standard

Tecna cod. SD933

|                      |   |
|----------------------|---|
| CATALOG #            | 300650  |
| SOURCE               | Isolated from <i>Microcystis aeruginosa</i>   |
| PRESENTATION         | Solution  |
| CONCENTRATION        | 10 ug/mL in Methanol, 1 mL  |
| PURITY               | >95% HPLC, 238 nm   |
| STORAGE              | Aliquot and store at < -20 °C   |
| STABILITY/SHELF LIFE | Two years if stored at < -20 °C   |
| APPEARANCE           | Clear   |
| HANDLING             | If you are not fully trained or are unaware of the hazards involved, do not use these compound. |

**FOR RESEARCH USE ONLY**