



Phone: (215) 357-3911

FAX: (215) 357-5232

WEB: Abraxiskits.com

## Desmethyl Microcystin-LR Standard

Tecna cod. SD931

CATALOG #	300656
SOURCE	Isolated from <i>Microcystis aeruginosa</i>
PRESENTATION	Solution
CONCENTRATION	10 ug/mL in Methanol, 1 mL
PURITY	>95% HPLC, 238 nm
STORAGE	Aliquot and store at < -20 °C
STABILITY/SHELF LIFE	Two years if stored at < -20 °C
APPEARANCE	Clear
HANDLING	If you are not fully trained or are unaware of the hazards involved, do not use these compound.

**FOR RESEARCH USE ONLY**