



Phone: (215) 357-3911
FAX: (215) 357-5232
WEB: Abraxiskits.com

Desmethyl Microcystin-RR Standard

Tecna cod. SD937

CATALOG #	300658
SOURCE	Isolated from <i>Microcystis aeruginosa</i>
PRESENTATION	Solution
CONCENTRATION	10 ug/mL in Methanol, 1 mL
PURITY	>95% HPLC, 238 nm
STORAGE	Aliquot and store at < -20 °C
STABILITY/SHELF LIFE	Two years if stored at < -20 °C
APPEARANCE	Clear
HANDLING	If you are not fully trained or are unaware of the hazards involved, do not use these compound.

FOR RESEARCH USE ONLY